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
The following pages were translated using Google Translate. For this reason, their correctness cannot be guaranteed. The German version of the operating instructions is legally binding.

Structural unit/work area:  
quantum optics

Activity: Handling a pressure vessel

### DESIGNATION

## Handling a pressure vessel

### DANGERS TO PEOPLE AND THE ENVIRONMENT

- Pressure vessels can be subjected to high pressures (up to 200 bar) and have massive destructive potential if they burst.
- Loose parts and loose fasteners can fly around when pressure increases.
- Leaks can result in the release of highly flammable or toxic substances or the creation of explosive gas mixtures in dangerous quantities.

### PROTECTIVE MEASURES AND RULES OF CONDUCT

- The procurement, storage, use and necessary protective measures for pressure vessels must be discussed in advance with the safety officer.
- Do not throw or impact pressure containers.
- Use pressure vessels only after instruction.
- Observe the operating instructions of the respective pressure vessel.
- Carry out a visual inspection for obvious defects and a leak test before use.
- Do not introduce explosive fuels into the pressure vessel.
- Make sure the seals are intact.
- Do not use grease or oil on the seals:
- Do not point your head or other body parts in the direction of a possible ejected part
- Closures of pressure vessels may only be opened once pressure has been equalized with the atmosphere.
- Pay attention to possible chemical effects of the pressure vessel contents.

### BEHAVIOR IN THE EVENT OF DISTURBANCES

- When releasing the gas, stop the pressure increase immediately.
- Release pressure.
- Correct the fault and press again.
- If the fault cannot be rectified, take the pressure vessel out of service immediately and arrange for repairs.
- Chemical effects due to pressure vessel contents → see chemical safety data sheet

### ACCIDENTS - FIRST AID - EMERGENCY CALL 112

- Make an emergency call.
- Initiate first aid measures, see below.
- Trained first aiders: see first aid notice

First aid services must be recorded in the first aid book.

Measures to be taken in case of toxic and/or aggressive gases or chemicals from the pressure vessel:

- After skin contact: Wash with plenty of water. Remove contaminated clothing. → Doctor
- After eye contact: Rinse opened eyelid for several minutes under running water, protecting the uninjured eye → ophthalmologist