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The following pages were translated using Google Translate. For this reason, their correctness cannot be guaranteed. The German version of the operating instructions is legally binding.

Structural unit/working area:
Quantum optics

Activity: Filling of liquid nitrogen and in-house transport/storage

DESIGNATION

Handling Dewar vessels or cryogenic containers

DANGERS FOR HUMANS AND THE ENVIRONMENT

Dangers to humans

- **Cold, frostbite and cold burns**

Body contact with liquid nitrogen can cause frostbite or cold burns (deep-cold liquefied nitrogen has a temperature of approximately -196 °C)

- **Increased risk of fire and explosion**

If air (liquid oxygen boils at -183 °C) comes into contact with liquid nitrogen, e.g. in open Dewar vessels, atmospheric oxygen can condense into the nitrogen. Over time, liquid oxygen can accumulate in the nitrogen. If this comes into contact with highly flammable materials, there is a risk of spontaneous ignition.

- **Risk of bursting due to pressure increase**

The inclusion of liquid nitrogen in containers not intended for this purpose without pressure equalization can cause the containers to burst at room temperature. Cold embrittlement occurs in plastic containers.

- **Lack of oxygen, risk of suffocation**

Particularly during filling and transfer processes, high nitrogen concentrations in the breathing air can occur unnoticed in rooms with poor ventilation due to the displacement of the breathing air. Possible acute risk of suffocation! (1 liter of liquid N_2 produces approx. 700 liters of gaseous N_2 when evaporated)

PROTECTIVE MEASURES AND RULES OF CONDUCT

Technical protective measures

- Equipment intended for direct handling of liquid nitrogen must be made of cold-resistant material.
- Rooms with filled Dewars must be adequately ventilated. Air inlets and outlets must not be closed. (Warning devices for oxygen deficiency are recommended).
- Dewar vessels that are filled with an overpressure of 0.5 bar require a safety valve - and/or a bursting disc. These vessels are subject to regular inspection (max. every 5 years) in accordance with DIN EN 1251-3.
- Have the inspection carried out/organized by VWR (O26/187).

Organisational protective measures

- Dewar vessels may only be filled with freezing media when they are clean and dry.
- Cryogenic containers (Dewars) must be transported in such a way that they cannot fall over.
- Non-pressurized containers containing liquid nitrogen may only be closed with a loosely fitting lid or plug so that pressure equalization is possible.
- These containers may only be filled and transported without pressure (loose lid!). (Dewars that are suitable for internal overpressure are marked accordingly). **Do not connect pressureless containers to extraction systems that require a tight connection.**
- Do not store liquid nitrogen in open containers for long periods of time
- When transporting liquid nitrogen in elevators, no persons are allowed to ride on board (chemical journey).

PROTECTIVE MEASURES AND RULES OF CONDUCT (continued)

Personal protective measures

- When filling and handling Dewar vessels, always use personal protective equipment, safety glasses and suitable protective gloves (material: Kevlar or special cryogenic gloves).
- If there is a risk of splashing, wear a face shield. Make sure that no liquid can get into shoes and gloves.

RESPONSE TO MALFUNCTIONS

Disturbances

- If large amounts of liquid nitrogen spill, leave the room immediately and warn other people.
- The room may only be entered with self-contained breathing apparatus until it can be proven that it is safe.
- Never carry out repairs yourself.
- Exposure to fire may cause container to rupture/explode.

MAINTENANCE AND REPAIR

Maintenance and repair

- Containers without test markings may no longer be filled.
- Tests are organized by VWR (please register vessels there)
- Repairs may only be carried out by authorized and trained personnel.

WHAT TO DO IN THE EVENT OF ACCIDENTS - FIRST AID - EMERGENCY CALL 112

- Make an emergency call.
- Initiate first aid measures.
- After skin contact, pour large amounts of lukewarm water over cold-burned areas of the body and then cover loosely with sterile dressing material.
- Remove unconscious persons to fresh air while protecting yourself. Keep warm and calm. If breathing stops, give artificial respiration.
- Trained first aiders: see first aid notice
- First aid services must be recorded in the first aid book.

CONSEQUENCES OF NON-COMPLIANCE

Injuries

- Improper handling can lead to serious injuries

Property damage

- Improper handling may result in containers bursting or fires.

Legal consequences

- Operating instructions are binding and represent written occupational safety - instructions for employees.
- Failure to comply may result in legal consequences.
- Failure to follow this instruction is a violation of the given instructions and will be punished accordingly.